Vol.3 No.2

RapidDetectionofMycobacteriumtuberculosisbyusingMPT64Immunochromatographic Assay

Dinesh Babu R*, Arun Kumar P, Heetandra Kumar Purohit and Shivashri C

Department of Microbiology, BKL Walawalkar Rural Medical College, Ratnagiri, Maharashtra, India

Abstract

Tuberculosis (TB) is one of the major infectious diseases in the developing countries and is more prevalent in India. The confirmatory identification of M. tuberculosis (MTB) is still being done by the conventional, laborious and time-consuming biochemical methods. Before this period of confirmation, the diseased patient is prone to encounter susceptible individuals; this clearly indicates the need for rapid, enhanced identification of Mycobacteria along with an effective drug sensitivity testing treatment. A specific antigen, MPT64 (major culture filtrate protein), differentiates the complex of M. tuberculosis from mycobacteria other than tuberculosis (MOTT) species. The present study evaluated the identification of Mycobacteria tuberculosis from MOTT by using SD Bioline kit and the findings were correlated with a conventional biochemical test and the sensitivity, specificity, and predictive values of this kit were also assessed. 80 (n=80) mycobacterial isolates from Sputum and BAL were tested by using MPT64 immunochromatographic assay and validated by standard microscopic and biochemical assays. In the assay there was no significant difference observed in identifying MTBC isolates between SD MPT64 TB Ag Kit® assay and Biochemical test. Thus, SD Ag MPT64 kit displayed 100% sensitivity, specificity, positive predictive and negative predictive values. Hence the kit is highly appropriate for use in TB diagnosis because of its simplicity, cost effectiveness and rapidity

Keywords: : Tuberculosis; Infectious diseases; Mycobacteria

Introduction Tuberculosis (TB) caused by the infection of bacillus Mycobacterium tuberculosis (MTB) characteristically affects the lungs (pulmonary TB) but also may affect other sites (extra pulmonary TB). It is spread by air and potentially fatal disease of human [1]. The people infected by HIV are highly prone to develop TB disease [1]. It has been estimated by the World Health Organization (WHO) in 2013 that there about 9 million (range, 8.6 million-9.4 million) cases globally. There are about 2.0-2.3 million incident cases of TB annually reported in India, this accounts for a greater number than any other country (i.e.) a fifth of new cases in the world [1]. There is also an alarming increase in number of drug-resistant cases of TB. Therefore, the need for rapid and correct identification of Mycobacteria and rapid drug sensitivity testing for effective treatment of the disease is increasing [2,3]. It is well known that though the automated culture systems Bactec 460, MBBacT (Bio Merieux, France), MGIT have reduced the time of culture,

they do not differentiate MTB and MOTT. The conventional biochemical methods are being used for the identification of M. tuberculosis (MTB). These tests are time consuming and raise the need for safety precautions; during this process the infected individual spread the infection to many susceptible individuals. The current objective of WHO is to find rapid identification and detection of drug resistance within 2 days by using employing nucleic acid amplification assays. This needs a specialized set up and trained personal, which may not be affordable for poor countries [4]. The Standard Diagnostics, Inc. (SD) (Yongin, Korea) developed an MPT64-based, simple and rapid immunochromatographic assay known as the SD Bioline TB Ag MPT64 RAPID® test (SD bioline kit). This lateral flow test has been shown to identify the M. tuberculosis complex from the MOTT, by impregnating the mouse monoclonal anti-MPT64 antibody on nitrocellulose membrane [1,5,6]. The specific antigen MPT64 (24 kDa) is one of the major culture filtrate protein encoded by RD2 in genes, which differentiates the M. tuberculosis complex from the mycobacteria other than tuberculosis (MOTT) species [7,8]. The actively secreted proteins of M. tuberculosis are the first to interact with the host immune system and recognized by human Th1 cells. Therefore, it could be useful for TB diagnosis or as part of a novel candidate vaccine against TB [9,10]. A large amount of variability in the diagnostic accuracy of MPT64 has been reported, depending on the recombinant antigen used in assays [7,8]. The present study used SD Bioline kit to evaluate the identification of MTB from MOTT and the findings were correlated with conventional biochemical test, and the sensitivity, specificity, and predictive values of this kit were also assessed.

Materials and Methods

80(n=80) mycobacterial isolates from Sputum and BAL, one control strain H3R7V, 15(n=15) Gram positive and gramnegative bacterial isolates from Urine, Sputum and BAL were tested by using MPT64 immunochromatographic assay and validated by standard microscopic and biochemical assays.

Immunochromatographic (ICT) assay

As per the manufacturer's instructions, MTB Antigen MPT64 rapid ICT kit (SD Bioline, Seoul, South Korea) was used [11]. The confirmation of MTB isolates was done by immobilizing the test kit with Mouse monoclonal antibodies, against MPT64 antigen, on a nitrocellulose membrane. The entire procedure

ThisworkispartlypresentedatJointEvent 9th Annual Congress onClinical Microbiology and Infectious Diseases

Vol.3 No.2

was done inside a biosafety class II cabinet. The 2-4 isolated colonies scraped from the solid LJ medium was suspended in 200 μ l of extraction buffer. Further this emulsified solution was applied in the sample well. After 15 minutes of incubation on sample application, the test was interpreted at room temperature. The entire procedure was carried out as per manufacturer instruction. The control region did not show any band, which was considered invalid. The test results were validated by biochemical test (Niacin test and Nitrate Reduction test).

Results

80(n=80) mycobacterial isolates were cultured in LJ media, after 8 weeks of incubation it was further confirmed by smear by using Ziel Neelson method. Validation of mycobacterial spp was done by niacin test and nitrate reduction test of these 80 isolates 75 were (n=75) MTB and 5(n=5) isolates were MOTT (Figure 1). All 80 isolates showed pink color band in control line and One H37Rv control strain shows pink colored band in both control and test line for confirming the presence of Tb MPT64 antigen. In the test line, 73 isolates of MTB displayed dark pink band, whereas 2 isolates showed faint band and were exposed to longer incubation period of 48 hours in the repeat test as well as biochemical test (niacin and nitrate reduction test) and was confirmed to be MTBC. It was clearly observed that there was no significant difference shown in detection of MTBC isolates between SD MPT64 TB Ag Kit® assay and Biochemical test. Hence, the sensitivity, specificity, positive predictive & negative predictive values of the SD Ag MPT64 kit was found to be 100%.

Conclusion

In this study, there was no significant difference found in the predictive values, sensitivity and specificity for isolates. The band intensity was more prominent in the liquid cultures. The study attempted to evaluate the ICT for the detection of MPT64 antigen in smear positive sputum samples. The total expenditure of ICT and biochemical tests will be approximately same. The swiftness of the test will markedly reduce the time of MTB culture and DST; thus, it can significantly contribute in TB control program. Our study results clearly demonstrated the use of SD Bioline test for the mycobacteria identification. The conventional identification tests can be effectively replaced by the immunochromatographic assays for detecting MP64 antigen. Thus, this rapid, simple and cost-effective method can be used in TB diagnostic laboratories.

Vol.3 No.2

ThisworkispartlypresentedatJointEvent 9th Annual Congress onClinical Microbiology and Infectious Diseases